Quill™ Self-Retaining System (SRS)
Composé de polypropylène Non-Absorbable
Surgical Suture Material.

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**Quill™ Self-Retaining System (SRS)
Comprised of Polypropylene
Non-Absorbable Surgical Suture Material**

DESCRIPTION

The Quill™ SRS comprised of polypropylene consists of barbed suture material, armed with a surgical needle on each end. Barbs allow for tissue approximation without the need to tie surgical knots.

Quill™ SRS comprised of polypropylene is an isotactic polypropylene polymer of high molecular weight. The material is provided dyed (blue) or undyed (clear). The pigment for the blue dyed suture material is Phthalocyaninato (2-) Copper.

While the formation of barbs in the Quill™ SRS comprised of polypropylene reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduce their effective strength. For this reason, the strength of the Quill™ SRS comprised of polypropylene can be compared to USP knot strength of non-barbed sutures.

TENSILE STRENGTH

USP SUTURE SIZE	SUTURE TENSILE STRENGTH (kgf)	QUILL™ SRS EQUIVALENT SIZE
1	2.72	2
0	2.16	1
2-0	1.44	0
3-0	0.96	2-0
4-0	0.60	3-0

ACTIONS

Quill™ SRS comprised of polypropylene elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the Quill™ SRS comprised of polypropylene by fibrous connective tissue. Quill™ SRS comprised of polypropylene is not absorbed, nor is any significant change in tensile strength retention known to occur *in vivo*.

INDICATIONS

Quill™ SRS comprised of polypropylene is indicated for use in soft tissue approximation excluding closure of the epidermis.

CONTRAINDICATIONS

Quill™ SRS comprised of polypropylene is not indicated for surface closures through the epidermis as the small opposing facing barbs make Quill™ SRS removal unfeasible.

WARNINGS

Do not resterilize. Discard opened, unused Quill™ SRS and associated surgical needles.

Users should be familiar with surgical procedure and techniques involving sutures before employing Quill™ SRS comprised of polypropylene, for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the *in vivo* performance (under ACTIONS section) when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.

The safety and effectiveness of Quill™ SRS comprised of polypropylene have not been established for use in fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissue, neural tissue, ophthalmic surgery, or for use in microsurgery, and the closure of other soft tissues that will be under excess tension.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

PRECAUTIONS

Quill™ SRS contains bidirectionally oriented barbs to anchor tissues and does not require knots to approximate opposing edges of a wound. Tying knots on the barbed section of the material will damage the barbs and potentially reduce the suture tensile strength and barb effectiveness. For the bidirectional forces to be created and for the device to function properly, both sides of the Quill™ SRS must be engaged in the tissue. Additionally, when completing placement, an additional backstitch or bite of tissue lateral to the end of the incision is required to lock the device in place.

Avoid contacting Quill™ SRS and associated needles with other materials (e.g. surgical gauze, drapes, etc.) in the surgical field to prevent ensnaring on the barbs. If the barbs catch, carefully pull the material in the opposite direction of the needle to disengage it from the barbs.

Care should be taken to avoid damage when handling. Avoid crushing or crimping the suture material with surgical instruments, such as needle holders and forceps. Do not pull Quill™ SRS out of the package by the needles as this can cause the barbs to catch on one another. Do not attempt to remove memory in the polymer by running fingers down the suture material as this can damage the barbs.

Infections, erythema, foreign body reactions, transient inflammatory reactions and in rare instances wound dehiscence are typical or foreseeable risks associated with any suture and hence are also potential complications associated with Quill™ SRS comprised of polypropylene.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of this device may include, wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, and pain, edema and erythema at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

APPLICATION

Use as required per surgical procedure. Deeper layers of the wound are to be closed in such a fashion as to take tension off of the dermal edges.

To deploy Quill™ SRS, the first end should be pulled through the tissue until the center transition zone has reached the tissue. Estimation of the center can be aided by taking a single bite of tissue, then aligning the two needles together. Taking one arm of the Quill™ SRS, at least two loose bites of tissue should be taken. This is then repeated with the other arm of the Quill™ SRS. Once at least four bites, or passes through the tissue have been completed, each of the strands can then be grasped and the tissue can be approximated to the desired tension.

To remove the entire Quill™ SRS comprised of polypropylene: Cut the Quill™ SRS at the midline, between the opposing barbed segments, and then pull the lateral ends to remove the device.

STERILITY

Quill™ SRS comprised of polypropylene is sterilized by ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged. Discard opened unused sutures. Do not use after expiration date.

HOW SUPPLIED

Quill™ SRS comprised of polypropylene is available sterile, in various barb configurations (e.g. 7 cm x 7 cm), dyed blue or undyed and USP diameter sizes 3-0 through 2 (metric sizes 2 through 5). Each Quill™ SRS comprised of polypropylene is double armed with various size needles.

SYMBOLS USED FOR LABELING

-  Product Code
-  Consult, see instructions for use
-  Do Not Reuse
-  Use by
-  Sterilized using ethylene oxide
-  Do Not Resterilize
-  Lot Code
-  Date of Manufacture
-  Manufacturer
-  Do not use if package is damaged

Rx Only **CAUTION:** Federal (USA) law restricts this device to sale and use by, or on the order of, a physician.

 CE mark and identification number of notified body. Product conforms to the essential requirements of The Medical Device Directive 93/42/EEC.

 Authorized Representative in the European Community.

PBN MEDICALS Denmark A/S
Knud Bro Alle 3, DK-3660 Stenløse
Denmark

To reorder product, please contact your local distributor, sales representative or Customer Service at:

Angiotech
100 Dennis Drive
Reading, PA 19606 USA
Telephone: 732 626-6466 • 877 991-1110 (USA and Canada Only) Fax: 610 404-4010

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QUILL™ is a trademark of Quill Medical, Inc. Quill™ is used
the trademark owner.
All Rights Reserved.