Quill™ Self-Retaining System (SRS)
Composé de Dyeed PDO (Polydioxanone)
Synthétique Absorbable Surgical Suture Material

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Rx Only CAUTION: Federal (USA) law restricts this device to sale and use by, or on the order of, a physician.

Quill™ is a trademark of Quill Medical, Inc. Quill™ is used by Surgical Specialties Corporation with permission from the trademark owner.

© 2009 Angiotech Pharmaceuticals, Inc.
© 2009 Surgical Specialties Corporation
All Rights Reserved.

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PBN MEDICALS Denmark A/S
Knud Bro Alle 3, DK-3660 Stenløse
Denmark

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**Quill™ Self-Retaining System (SRS)
Comprised of Dyed PDO (Polydioxanone)
Synthetic Absorbable Surgical Suture Material**

DESCRIPTION

The Quill™ SRS comprised of dyed PDO consists of barbed suture material, armed with a surgical needle on each end. Barbs allow for tissue approximation without the need to tie surgical knots.

Quill™ SRS is comprised of dyed (violet) polyester, poly(p-dioxanone). The empirical molecular formula of which is (C₄H₆O₃)_x. The pigment for the violet dye is D&C Violet No. 2. Polydioxanone has been found to be nonantigenic, nonpyrogenic and to elicit only a slight tissue reaction during absorption.

While the formation of barbs in the Quill™ SRS comprised of dyed PDO reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduce their effective strength. For this reason, the strength of the Quill™ SRS comprised of dyed PDO can be compared to USP knot strength of non-barbed sutures. Additionally, USP designations for diameter are used to describe the Quill™ SRS comprised of dyed PDO prior to barbing, except for minor variation in suture diameter with a maximum overage of 0.1 mm.

TENSILE STRENGTH

USP SUTURE SIZE	SUTURE TENSILE STRENGTH (kgf)	QUILL™ SRS EQUIVALENT SIZE
5-0	0.68	4-0
4-0	0.95	3-0
3-0	1.77	2-0
2-0	2.68	0
0	3.90	1
1	5.08	2

ACTIONS

Two important characteristics describe the in vivo performance of absorbable sutures: first, tensile strength retention and second, the absorption rate (loss of mass). Quill™ SRS comprised of dyed PDO has been formulated to minimize the variability of these characteristics and to provide wound support through the critical wound healing period and an extended healing period.

The results of implantation studies in animals using Quill™ SRS comprised of dyed PDO indicate that for sizes 3-0 and larger, approximately 80% of the original strength remains after two weeks and four weeks of implantation. At six weeks post implantation, approximately 40% to 70% of the original strength is retained. For size 4-0, approximately 67% of the original tensile strength remains after two weeks, approximately 50% at 4 weeks and approximately 37% at 6 weeks.

Data obtained from implantation studies show that absorption of Quill™ SRS comprised of dyed PDO is minimal until about 120 days, and is essentially complete within 180 days (six months).

INDICATIONS

Quill™ SRS comprised of dyed PDO is indicated for use in soft tissue approximation where use of absorbable sutures is appropriate.

CONTRAINDICATIONS

Quill™ SRS comprised of dyed PDO is not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and is not to be used in conjunction with or for fixation of prosthetic devices (e.g. heart valves or synthetic grafts) that are non-absorbable in nature.

WARNINGS

Do not resterilize. Discard opened, unused Quill™ SRS and associated surgical needles.

Users should be familiar with surgical procedure and techniques involving absorbable sutures before employing Quill™ SRS comprised of PDO for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under ACTIONS section) when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.

The safety and effectiveness of Quill™ SRS comprised of dyed PDO has not been established for use in fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissue, neural tissue, ophthalmic surgery or for use in microsurgery, and the closure of other soft tissues that will be under excess tension. Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the surgeon.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, Quill™ SRS comprised of dyed PDO may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention or which may require additional support.

PRECAUTIONS

Quill™ SRS contains bidirectionally oriented barbs to anchor tissues and does not require knots to approximate opposing edges of a wound. Tying knots on the barbed section of the material will damage the barbs and potentially reduce the suture tensile strength and barb effectiveness. For the bidirectional forces to be created and for the device to function properly, both sides of the Quill™ SRS must be engaged in the tissue. Additionally, when completing placement, an additional backstitch or bite of tissue lateral to the end of the incision is required to lock the device in place.

Avoid contacting the Quill™ SRS and associated needles with other materials (e.g. surgical gauze, drapes, etc.) in the surgical field to prevent ensnaring on the barbs. If the barbs catch, carefully pull the material in the opposite direction of the needle to disengage it from the barbs.

Care should be taken to avoid damage when handling. Avoid crushing or crimping the suture material with surgical instruments, such as needle holders and forceps. Do not pull the Quill™ SRS comprised of dyed PDO out of the package by the needles as this can cause the barbs to catch on one another. Do not attempt to remove memory in the polymer by running fingers down the suture material as this can damage the barbs.

Infections, erythema, foreign body reactions, transient inflammatory reactions and in rare instances wound dehiscence are typical or foreseeable risks associated with any suture and hence are also potential complications associated with Quill™ SRS comprised of dyed PDO.

When using Quill™ SRS comprised of dyed PDO subcutaneously, the device should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

ADVERSE REACTIONS

Adverse effects associated with the use of this device may include, wound dehiscence, failure to provide adequate wound support in closure of the site where expansion, stretching or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from condition which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs and transitory local infection as the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

Due to prolonged suture absorption, some irritation and bleeding may occur.

APPLICATION

Use as required per surgical procedure. Deeper layers of the wound are to be closed in such a fashion as to take tension off of the dermal edges.

To deploy Quill™ SRS, the first end should be pulled through the tissue until the center transition zone has reached the tissue. Estimation of the center can be aided by taking a single bite of tissue, then aligning the two needles together. Taking one arm of the Quill™ SRS, at least two loose bites of tissue should be taken. This is then repeated with the other arm of the Quill™ SRS. Once at least four bites, or passes through the tissue have been completed, each of the strands can then be grasped and the tissue can be approximated to the desired tension. Bites, or passes through the tissue, can be taken to approximate the wound. Care should be taken to utilize Quill™ SRS on the barbed segments only. Do not attempt to approximate wounds using the non-barbed segments, as the barbs are required for successful wound approximation with Quill™ SRS.

For Quill™ SRS with straight needles, Quill™ SRS can be deployed at right angles to the wound. Starting at one edge of the wound, the Keith needle is deployed at right angles to the wound along the desired tissue plane and the needle is advanced through the tissue until the desired amount of tissue has been engaged, then the needle can exit at the desired location remote from the wound edge. The same procedure is repeated with the second arm of the Quill™ SRS. After successful deployment of both arms of the Quill™ SRS and needles exits, the center transition zone of the Quill™ SRS should be roughly in the center of the open wound. The ends of the Quill™ SRS where they exit remote to the wound can now be grasped and the tissue can be pushed together to facilitate the soft tissue approximation to the desired tension.

To remove the entire Quill™ SRS comprised of dyed PDO: Cut the Quill™ SRS at the midline, between the opposing barbed segments, and then pull the lateral ends to remove the device.

STERILITY

Quill™ SRS comprised of dyed PDO is sterilized by ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged. Discard opened unused sutures. Do not use after expiration date.

HOW SUPPLIED

Quill™ SRS comprised of dyed PDO is available sterile, in various barb configurations (e.g. 7 cm x 7 cm), dyed violet in various lengths and USP diameter sizes, 4-0, through 2 (metric sizes 1.5 through 5). Quill™ SRS comprised of dyed PDO is double armed with various size needles.

SYMBOLS USED FOR LABELING

- Product Code
- Consult, see instructions for use
- Do Not Reuse
- Use by
- Sterilized using ethylene oxide
- Do Not Resterilize
- Lot Code
- Date of Manufacture
- Manufacturer
- Do not use if package is damaged
- CAUTION:** Federal (USA) law restricts this device to sale and use by, or on the order of, a physician.
- CE mark and identification number of notified body. Product conforms to the essential requirements of The Medical Device Directive 93/42/EEC.
- Authorized Representative in the European Community.
PBN MEDICALS Denmark A/S
Knud Bro Alle 3, DK-3660 Stenløse
Denmark

To reorder product, please contact your local distributor, sales representative or Customer Service at:

Angiotech
100 Dennis Drive
Reading, PA 19606 USA
Telephone: 732 626-6466 • 877 991-1110 (USA and Canada Only) Fax: 610 404-4010