**P**

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# UltraFit™ PK System

## Coronet Donor Punch Set

Surgical Specialties Corporation

100 Dennis Drive

Reading, PA 19606 U.S.A.

Phone: 610-404-1000

www.angiotech.com

Manufactured For:

Surgical Specialties Corporation

100 Dennis Drive

Reading, PA 19606 U.S.A.

Phone: 610-404-1000

www.angiotech.com

SS0024.2

## USA SYMBOLS

<b>LOT</b> Lot Number	 Use Until Date
<b>STERILE EO</b> Sterilized by Ethylene Oxide. Do not use if product is opened or damaged.	 Do Not Re-use

<b>REF</b> Catalog Number	 See Insert
---------------------------	--

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## USA

### Coronet Donor Punch Set Instructions for use

#### Description.

The Coronet Donor Punch Set is a disposable punch and trephine system, for the trephination of the donor cornea button.

#### Intended Use.

The Coronet Donor Punch Set allows the surgeon to mark the four cardinal suture points and holds the donor corneal tissue in the most relaxed state in preparation for the cutting of the donor button. The trephine, supplied as part of the system, is ultra sharp and makes a clean vertical cut through the corneal tissue.

#### CAUTION

- This device is intended for SINGLE USE ONLY. It is not recommended for repeat use.
- The punch is supplied sterile and ready for use. Do not use if the pouch has been opened or damaged.
- The Coronet trephine is much sharper than a traditional trephine. Caution should be exercised at all times when handling and using the trephine.
- The Coronet trephine is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard the blade if damaged.

#### Sterilization

- The trephine is supplied sterile and ready for use.
- DO NOT RE-STERILIZE.**

#### Surgical Technique

- Separate the two parts of the punch – the base and the trephine guide.
- A gentian ink marking pen can be used to mark the guide holes for the cardinal suture points.
- Remove the cornea from the storage medium by grasping the rim of the surrounding sclera with toothed forceps and dry the anterior surface with a lint free PVA Sponge spear.
- Position the cornea in the block well with the endothelial surface up. Center the cornea using the guide holes.
- To replace the trephine guide, position the guide over the base so that it will drop vertically and press down gently until the two parts are correctly located. Make certain that the trephine guide is pressed down fully onto the base before proceeding. If the two parts do not locate correctly, either rotate the trephine guide counter-clockwise, lift it off the base and reposition or press down firmly on the top of te guide with the thumb or finger until the two parts align correctly.
- Remove the trephine from the packaging and position in the guide, taking care not to damage the ultra sharp cutting edge.
- Confirm the position of the corneal tissue and if necessary use a micro hook instrument via the windows of the trephine guide, to manipulate and position the corneal tissue optimally.
- Lower the trephine blade gently onto the endothelial surface, and then press with a gentle but firm movement on the trephine to achieve the vertical cut through the corneal tissue. It is not necessary to rotate the trephine but if in doubt the trephine and coronet can be rotated together to ensure a thorough cut.
- Please note. The exceptionally super sharp trephine blade will cut through the corneal tissue without any sound or resistance. Once resistance is felt, the cut has been completed.**
- Once the cut has been achieved, remove the trephine and the trephine guide together in one piece. The corneal button, peripheral cornea and scleral ring will remain in the well of the base.
- In the unlikely event that the corneal button remains in the lumen of the blade, a balanced salt solution or a viscoelastic substance should be gently dripped into the center of the blade to dislodge the button.

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